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
Key appearance	Key purpose
	The key is used to turn on and off the pH-meter. Holding time for response – 2 s.
	a) Measurement mode: to turn to the pH-meter calibration mode. Holding time for response – 0.5 s. b) Work with the scratchpad and display menus: to move leftward in a line.
	a) Measurement mode: to select mode of pH, pH ₂₅ , EMF measurement (out of available modes, see cl. 1.5.8). Holding time for response – 0.5 s. b) Work with the scratchpad and display menus: to move upward in a line.
	a) Measurement mode: to turn from measurement mode to the mode of viewing records in the scratchpad. Holding time for response – 0.5 s. b) Work with the scratchpad and display menus: to move rightward in a line.
	a) Measurement mode: to turn on and off display light. b) Work with the scratchpad and display menus: to exit from the scratchpad and display menus.
	a) Measurement mode: to enter the display menu. Holding time for response – 0.5 s; b) Work with the scratchpad and display menus: to move downward in a line.
	a) Measurement mode: to enter data into the scratchpad. Holding time for response – 0.5 s. b) Work with the scratchpad and display menus: to confirm preset parameters and selected operating modes.

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